

## Using the Manage XchangER Configurations

The Manage XchangER Configurations page will allow an XchangER Admin to manage configuration options for XchangER. These options include:

1. Install an authentication cookie or authorize a device's public IP address to provide designated devices access to the XchangER runs panel.
2. Modify the Facility Name which is displayed at the top of the runs panel.
3. Include ePCR's on your runs panel when your facility has been designated a Transport To Facility (on by default), an On-Line Medical Facility, or a Transfer To facility.
4. Include the Patient Initials in the Runs Panel.
5. Include the Booklet Number in the Runs Panel and Filter criteria.
6. Set the default number of PCR's per page.

To begin the authentication process on a new device, navigate to <https://xchanger.healthems.com>. A login window will be displayed. Users with the XchangER Admin role will be able to use their credentials to access the Install XchangER Configurations page.



**XchangER™**  
for HEALTHEMS®

**Login**  
Please log in.

Username:

Password:



[Did you forget your Password?](#)

On a previously authenticated device which has access to the XchangER runs panel, click the wrench icon in the top right corner to be prompted to log in as an XchangER Admin to access the Install XchangER Configurations page.

Print	Ho	Svc Date	Left	Si	Source	Cc	Call #	Agency	Veh	PT	Att	Age	M/F	Dispatch	Descript	Chief Complaint	Provider Impressions	Vitals	Initial Assessment	
		08/10/2011	07:51	ePCR			872223	HealthEMS Ambu	E11											
		08/10/2011	07:51	ePCR			872222	HealthEMS Ambu	E11		BDJ					my fingers hurt				
		08/08/2011	16:29	ePCR			101007	HealthEMS Ambu	M18	JDD		2y	M		Chest Pain	Chest Pain	Cardiac/Chest Pain: Chest Pain Cardiac Symptoms Cardiac/Chest Pain: Cardiac Sym Chest Pain	Time=16:14; BS=88; BP=120/80; CO2=35; Pulse=68; Resp=16; SPO2=96;	GC8=15; GC9=15; AVPU=Alert; AIR=Patent BREATH=Normal-Unlabored; RHYTHM=NSR; RHYTHM2=NSR;	
		07/25/2011	10:22	ePCR			555555	HealthEMS Ambu	1094							my fingers hurt				
		05/27/2011	11:05	ePCR			6002	HealthEMS Ambu	M18	JD		65y	M		Chest Pain	Chest Pain	Hyperglycemia Pancreatitis: Acute Cardiac/Chest Pain: Chest Pain Possible Intoxication-Check the MOI Abnormal Involuntary movements Cardiac Symptoms Cardiac Arrest Cardiac/Chest Pain: Cardiac Sym Back pain, traumatic, lumbar or tho Embolism: Pulmonary Urinary Retention Asthma Attack Cardiac/Respiratory Arrest: Card Confusion Airway Obstruction Chest Pain PAINFUL RESPIRATION	Time=11:02; BS=96; BP=116/82; Pulse=96; Resp=12; SPO2=96;	GC8=15; GC9=15; AVPU=Alert; AIR=Patent BREATH=Normal-Unlabored; RHYTHM=NSR; RHYTHM2=NSR;	

On the Install XchangER Configurations page, click the Configure button to go to the Manage XchangER Configurations page.

Available Configurations	Config Id	Description	Install Cookie	Configure
	123456	XER - St. Mary's Medical Center		

On the Manage XchangER Configurations page, an XchangER Admin may add the public IP address of a pc by clicking the “Add Ip Address” link. XchangER will automatically populate the current public IP address being displayed by the pc. As an additional security measure, if the check box for “Restrict access to configured IPs” is checked, the Install Cookie button will only appear if the pc is displaying an approved public IP address.

XchangER Client Configuration

Facility Management

Facility Name: St. Mary's Medical Center

Show Booklet Number

Include Transfer To Runs

# Results Per Page (Default): 10

Show Patient Initials

Include On-Line Medical Facility Runs

Include Transport To Runs

System Access

Restrict access to configured IPs  **Install Cookie**

[Add Ip Address](#)

**Submit**

[Navigate to XchangER Runs Page](#)

To save any changes made to the Manage XchangER Configurations page, click Submit. Users will be required to enter an XchangER Admin credential to save the changes.

Please log in to update the configurations.

Username:

Password:

Did you forget your Password?

**Log In** **Cancel**

When finished, click the “Navigate to [XchangER Runs Page](#)” link at the bottom of the Manage XchangER Configurations page to return to the main XchangER runs panel.